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USER MANUAL FOR THE RESULTS OF GENETIC  
PREDISPOSITION TEST TO DEVELOPING AN HEPATO-  
CELLULAR CARCINOMA

Results generated by the in vitro diagnostic medical device  
(software) manufactured by :

**BILHI**  
G E N E T I C S

BILHI GENETICS, 60 Avenue André Roussin - 13321 MARSEILLE  
Cedex 16  
e-mail : support@bilhigenetics.com

In partnership with the biomedical laboratory Gen-Bio  
Gravanches - 8 rue Jacqueline Auriol - 63100 Clermont-Ferrand  
/ Tel. : +33(0)4 73 23 25 84

#### PRODUCT

BILHI LIVER HCC SOFTWARE is an in vitro diagnostic medical device, which determines the genetic predisposition to developing an hepato-cellular carcinoma (HCC), based on genetic data obtained from human biological sampling. It has been designed and developed by BILHI GENETICS as a software.

#### DESCRIPTION

This software uses genetic data to generate a numerical statistical result in order to evaluate the genetic predisposition to developing an hepato-cellular carcinoma.

#### INDICATIONS, CONTRAINDICATIONS AND POTENTIAL UNDESIRE EFFECTS :

##### Indications

The use of the software results is indicated for patients with hepatic cirrhosis caused by the hepatitis C virus and willing to know the risk of developing an hepato-cellular carcinoma. Patients should be caucasian and over 20 years old.

**WARNING: THIS TEST MUST ONLY BE USED IF A MEDICAL PRESCRIPTION HAS BEEN PROVIDED BY THE PRACTITIONER AND INFORMED CONSENT HAS BEEN GIVEN BY THE PATIENT**

WARNING: Under French Law, the sale of genetic tests in France is restricted to medical practitioners and their use for preventative medicine or scientific research.

##### Contraindications

The (non-exhaustive) contraindications are the following:

1. Mental illness
2. Patient under 20

3. Any other medical complication that would compromise a following successful sample and/or the results of the device
  4. Patient of non-Caucasian ethnicity
- Software results are not designed, intended or sold for any other purpose than those listed above.

##### Potential side effects

A list of undesirable effects related to the use of the software results, although not exhaustive, is the following :

2. Behavioural change in the patient regarding his/her decision whether or not to undergo treatments proposed by practitioner

PLEASE NOTE: Depending on the results of the software, a therapeutic adaptation can be decided upon by the practitioner who receives the final exam report.

##### USE OF SOFTWARE RESULTS

**WARNING: The results generated by the software are neither for diagnostic purposes, nor should be used to predict the appearance of a pathology or symptoms, nor to detect or identify any specific treatment.**

**The results are given in the form of a patient classification in one of seven distinct following groups:**

- **Group 0.1** meaning a **lower risk** than in the general population risk group: **2%** of patients in this group could develop an HCC

- **Group 1.0** meaning an **average risk** as in the general population risk group: **25%** of patients in this group could develop an HCC

- **Group 1.6** meaning a **1.6 times increased risk** than in the general population risk group: **40%** of patients in this group could develop an HCC

- **Group 2.4** meaning a **2.4 times increased risk** than in the general population risk group: **59%** of patients in this group could develop an HCC

- **Group 2.7** meaning **2.7 times increased risk** than in the general population risk group: **67%** of patients in this group could develop an HCC

- **Group 3.2** meaning **3.2 times increased risk** than in the general population risk group: **80%** of patients in this group could develop an HCC

- **Group 4.0** meaning **4 times increased risk** than in the general population risk group: **100%** of patients in this group could develop an HCC

**WARNING: A result deeming the patient to be at an increased risk does not necessarily mean that he/she will develop an HCC. A result deeming the patient to be at a decreased risk does not mean that he/she will not develop an HCC.**

##### WARNINGS AND PRECAUTIONS

Every evaluation made by the software will not necessarily lead to a result. It is possible that the software is not able to determine the genetic predisposition, especially if the procedures for collecting and preparing the genetic data have

not been respected. In this case, a new sample will be taken free of charge to confirm the previous result. Compliance with sampling procedures, including sample submission and identification as well as the preparation of genetic data files, are important factors towards the success of the software. In addition, patient cooperation during the sampling will greatly influence the results. If sample is saliva based, patients who smoke prior to specimen collection have been shown to be at risk of DNA destruction.

##### OTHER PREANALYTIC AND POSTANALYTIC WARNINGS ARE AS FOLLOWS

##### Selection of the genetic test

It is important to choose to prescribe genetic test BILHI LIVER HCC for patients with hepatic cirrhosis caused by the hepatitis C virus for an early management of a possible development of an HCC.

##### Preanalytic

Only patients who meet the criteria described in the "Indications" should be selected. Patients meeting the criteria described in the "Contraindications" should not be selected.

##### Postanalytic

The physician must provide the patient with detailed instructions on the use and limitations of the software results. The patient should be given advice and explanation about his risk of developing an HCC.

Any results returned in error must be returned immediately to the biomedical laboratory where the examination report was transmitted from.

##### ADDITIONAL INFORMATION

For further information on this test, or any additional technical information, please contact BILHI GENETICS Customer Service or the biomedical laboratory. Contact details have been provided on this notice.

##### PRODUCT RECLAMATION

Any health professional (for example, any client or healthcare system user) with a claim or who is unsatisfied with the quality of the software's results, identity, reliability, safety, efficiency and/or performance, should notify the biomedical laboratory or BILHI GENETICS.

For any claim, please provide the name of the genetic test, the patient identification number, as well as your name, address and the nature of your claim.

##### INFORMATION FOR THE PATIENT

The BILHI LIVER HCC genetic test is a test for a genetic predisposition to an hepato-cellular carcinoma as well as a preliminary and preventive measure for an early management by the practitioner. It is carried out via a partnership between the biomedical laboratory and BILHI GENETICS, who combine their respective technology, using your DNA sample upon your consent.

The BILHI LIVER HCC SOFTWARE uses your genetic data transmitted by the partner medical analysis laboratory. This software is uniquely intended for the BILHI LIVER HCC genetic predisposition test. The specific indications of the BILHI LIVER HCC SOFTWARE are the following :

- Assessment of the risk of progression of hepatic cirrhosis caused by the hepatitis C virus to hepatocellular carcinoma.

Early management of a risk of hepato-cellular carcinoma may lead to additional treatments to reduce the established risk. Therefore, it is very important that you carefully follow the advice of your medical practitioner upon reception of the results.

Following these instructions, you will increase your chances of getting the best care management.

##### WITHDRAWAL AND DISPOSAL OF GENETIC DATA

The removal and disposal of the genetic data once it has been treated by the software take place 30 years following the end of the software's "life", unless the patient withdraws his consent before that time by claiming his rights according to the General Regulations on the Protection of Personal Data. The removal and disposal of the genetic data takes place according to BILHI GENETICS internal procedure. No specific arrangements are necessary.

##### WARNING

The responsibility of the software manufacturer is limited to the applications and uses mentioned in this instruction manual. The use of results for any scientific or investigational use is strictly forbidden.

The results of the test for a genetic predisposition to hepato-cellular carcinoma have been generated by an In Vitro Diagnostic Medical Device

This In Vitro Medical Device has been evaluated according to a European directive 98/79. The results were deemed compliant. CE mark obtained : 03/2021

You must consult the instruction manual before using the results provided by the software.